BANiTi E35 BANiTi E45

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Operating Instructions

Table of contents

1	Before you begin	4
1.1	Structure of the document	4
1.2	Service life	7
2	Safety information	8
3	Technical description	10
3.1	Task.....	10
3.2	Contra-angle handpiece structure	10
3.3	Technical data	12
4	Preparation	15
4.1	Initial start-up and longer breaks in use	15
5	Operation	16
5.1	Replacing the instrument.....	16

5.2	Attaching and removing nickel titanium files	17
5.3	Setting the torque	18
5.4	Attaching the APEX clip	19
5.5	Connecting the instrument to an endometer	20
5.6	Recommendation for the treatment procedure.....	21
6	Re-processing.....	23
6.1	After each treatment session.....	23
6.2	At the end of the work day.....	23
7	Care and maintenance	24
7.1	Cleaning and maintaining the mechanics using spray	24
7.2	Cleaning the pushbutton chuck	26
7.3	Manually cleaning and disinfecting the surface	27
7.4	Sterilizing	29

Intended use

1 Before you begin ...

BANiTi E35 / BANiTi E45 is used for rotary processing of tooth substance for endodontics.

BANiTi E35 / BANiTi E45 complies with the latest regulations reflecting state-of-the-art technology. BANiTi E35 / BANiTi E45 meets standard ISO 7785 - 2.

1. Read the operating instructions prior to using BANiTi E35 / BANiTi E45.
2. Use BANiTi E35 / BANiTi E45 only for applications described in standard ISO 7785 - 2 and the operating instructions.
3. Observe the applicable hygiene standards, occupational safety regulations and accident prevention measures for BANiTi E35 / BANiTi E45.

Target group

This product is intended only for use by trained dental personnel.

1.1 Structure of the document

1.1.1 Identification of warnings

Warnings

- To prevent any personal injury, please observe all warnings.

Warnings can be identified as follows:

Instructions for use

 **DANGER!** indicates a danger leading to death or serious injury if not avoided.

 **WARNING!** indicates a danger that may lead to death or serious injury if not avoided.

 **CAUTION!** indicates a danger that may lead to injury if not avoided.

➤ To prevent material damage and additional expenses, please observe all instructions for use.

Instructions for use can be identified as follows:

 **NOTICE!** indicates measures for the prevention of material damage.

IMPORTANT: indicates information on the avoidance of additional expenses and other important information.

Tip: indicates information on making work easier.

1.1.2 Formats and characters

The symbols and character formats used in the present manual have the following meaning:

<ul style="list-style-type: none"> ✓ Prerequisite 1. First action step 2. Second action step or ➤ Alternative action ↪ Result 	Requests you to do something.
Use of formats and characters [→ 6].	Identifies a reference to another text passage and indicates the relevant page number.
• List	Identifies a list.

1.2 Service life

When used as intended:

- Non-moving parts have a typical service life of 5 years
- Moving parts have a typical service life of 3 years

No warranty claim can be inferred here, as wear may occur earlier or later than indicated above depending on use, frequency of sterilization and frequency of maintenance.

Elastomers, e.g. O-rings, must be replaced depending on their degree of wear.

Check the tool holding system annually for safety reasons.

Obligations of the user

Preventing the spread of infection

Preventing cross-contamination

Magnetic fields

2 Safety information

- Use only equipment that is in full working order.
- Protect yourself, patients and others against any foreseeable dangers. In order to do so, follow the safety instructions.
- Use the equipment as intended.
- You should always keep these operating instructions within reach for further reference.

Sterilize the instrument after each patient. This prevents any transmission of infective agents such as HIV (the AIDS virus) or hepatitis B which might cause serious illnesses or even lead to death.

Prevent cross-contamination between patients, users and third parties as follows: Sterilize the instrument after each patient.

Take the appropriate hygiene measures e.g. wear protective gloves.

The BANiTi E35 / BANiTi E45 instrument generates magnetic fields which might affect cardiac pacemakers. Therefore, the BANiTi E35 / BANiTi E45 instrument must not be used on patients or operated by users who have cardiac pacemakers!

Direct contact of the instrument with magnetic storage media can lead to data loss.

The instrument is therefore marked with the following symbol:

Instrument head overheating

If the instrument is defective, the area around the instrument head may heat up, thus creating a risk of burning the patient's oral mucosa.

Malfunction or damage

Immediately discontinue use in the case of malfunction or damage. Damaged instruments may cause injury. Notify the dental depot or the manufacturer.

Repair

Do not repair the instrument yourself.

The BA instrument is the result of a precise development and must not be disassembled by third parties. In the case of service or repair work, please send the instrument in an assembled state to B.A. International Ltd. or to a service center approved by B.A. International.

Accessory parts

Use only original accessory parts produced by the manufacturer.


If you have any questions, please contact your dental depot or the manufacturer.

3 Technical description


3.1 Task

The instrument serves for transmitting the driving power and speed of the electric motor or air motor (handpiece coupling according to ISO 3964) to the preparation tool.

3.2 Contra-angle handpiece structure

A	Pushbutton
B	Opening of chuck system
C	Gear ratio
D	Torque levels
E	Control ring
F	Joint for clip

Accessories

G	Insulating sleeve
H	NiTi file
I	Clip Ø 2.0 mm
J	Clip Ø 1.5 mm

Torque CARD

For the file assortment of each manufacturer, a corresponding Torque CARD is included in the accessories.

Here you can read the torque levels (K) recommended for the files used.

H658 M 675LN 1 / \$ S 07 00001

(A)
(B)
(C)

HSI HIBC code

A	Product label code (in this case: 675LN 1)
B	Year of manufacture (in this case: 2007)
C	Serial number (in this case: 00001)

3.3 Technical data

Instrument

	BANiTi E35	BANiTi E45
Gear ratio	115:1	66:1
Maximum motor speed in rpm	40000	40000

		BANiTi E35	BANiTi E45
Maximum operating speed in rpm (must be chosen according to the file manufacturer)		approx. 350	approx. 600
Motion		360° (fully rotating)	360° (fully rotating)
Pushbutton cover clamping system		NiTi files	NiTi files
Clips available for standard endometer connections			
	ø 1.5mm	x	x
	ø 2mm	x	x

Nickel titanium files

	BANiTi E35 / BANiTi E45
Shank standard	ISO 1797-1
Shank diameter in mm	2,334 - 2,350
Maximum total length in mm	25
Maximum working diameter in mm (ISO 021)	2,1

Coupling

	BANiTi E35 / BANiTi E45
Standard	ISO 3964

4 Preparation

4.1 Initial start-up and longer breaks in use

- Sterilize the instrument and accessories prior to startup.

5 Operation

⚠ CAUTION! Prior to the treatment, ensure that the file is securely positioned. A loose or partially removed file can detach itself from the head or break off. This leads to a risk of injury! Therefore, only use the instrument if the file is at least 10 mm deep and clamped securely in place.


⚠ CAUTION! Do not pull the patient's cheek back with the contra-angle handpiece while the motor is running. This would actuate the pushbutton, thus creating a risk of burning the patient's oral mucosa.

5.1 Replacing the instrument

⚠ CAUTION! The instrument should only be fitted or removed when the motor is at standstill.

Attaching the BANiTi E35 / BANiTi E45

- ✓ The motor has come to a stop.
- Insert the instrument until it snaps into place.

5.2 Attaching and removing nickel titanium files

⚠ NOTICE! Use only faultless NiTi files in order to prevent fatigue fractures of the files as far as possible.

Please take the following points into consideration when selecting files:


- The NiTi file is designed for fully rotating use.
- A torque CARD is available for the NiTi file.

Attaching a file

- ✓ The drive has come to a stop.
1. Insert the file by gently rotating it until it snaps into place. In doing so, do not press the pushbutton.
 2. Then pull and turn the file to check that it is firmly seated.


Removing a file


- ✓ The file has come to a stop.
- Remove it while pressing and holding down the pushbutton.

5.3 Setting the torque

- ✓ The required NiTi file is inserted [→ 17].
 - ✓ The speed set at the supply unit corresponds to the file manufacturer's specifications.
1. Read out the recommended torque level from the torque CARD of the file used.
 2. Turn the ring (A) until its front edge is located opposite the mark corresponding to the desired level 1 - 5.


5.4 Attaching the APEX clip

for connecting e.g.:

- Apex Pointer (Micro Mega)
- iPex (NSK)
- Mini Apex Locator (SybronEndo)
- ProPex (Dentsply Maillefer)
- ProPex II (Dentsply Maillefer)
- RayPex 4 (VDW)
- RayPex 5 (VDW)
- Root ZX (Morita)

for connecting e.g.

- Apex Master (EMS)

5.5 Connecting the instrument to an endometer

The BANiTi E35 / BANiTi E45 makes it possible to connect a measuring cable for an endometer by means of a clip in order to check the exact position of the file in the root canal.

⚠ NOTICE! The endometers must comply with the requirements of IEC 60601.

1. Please observe the operating instructions pertaining to the endometer being used.
2. Attach an APEX clip [→ 19].
3. Connect the measuring cable of the endometer to the clip's contact pin.
4. Perform the measurement in accordance with the equipment manufacturer's specifications.

Based on in vitro studies, we recommend defining the maximum preparation depth at the coronal end of the apex interval on the display of the endometer used.

Calibration


- ✓ The manufacturer of the endometer stipulates that a calibration must be performed.


1. Insert a file into the contra-angle handpiece [→ 17].
2. Connect the cable of the file clamp (without the actual file clamp) onto the clip of the contra-angle handpiece.
3. Perform the calibration in accordance with the equipment manufacturer's specifications.

5.6 Recommendation for the treatment procedure

⚠NOTICE! Use only faultless NiTi files in order to prevent fatigue fractures of the files as far as possible.

- The different file systems require different working techniques. For this reason, please always observe the information provided by the file manufacturer.
- Work without exerting any pressure on the file.
- Prior to the first treatment, we recommend performing conditioning tests on resin blocks or extracted teeth.
- Due to the special characteristics of the nickel-titanium material, the crown down technique should be preferred when preparing the root canal.

- When performing the APEX measurement, wear appropriate insulated gloves. We recommend performing the treatment with a cofferdam, while using an insulating sleeve. This will prevent inaccurate measurements caused by undesired leakage currents. The root canal file must not contact the patient's mucous membrane or the mucosal electrode during the measurement.
- If the limit torque is exceeded, a gentle creaking may be heard and vibrations felt. Remove the file by gently pulling it toward the coronal end. If this is not possible, the drive should be switched to counterclockwise rotation and the file withdrawn from the channel.
IMPORTANT: Subsequently, the file must be checked for damage or deformation (untwisting) and then replaced, if necessary.

6 Re-processing

6.1 After each treatment session

1. Clean and disinfect the instrument and accessories.
2. Apply spray to the instrument.
3. Sterilize the instrument and accessories.

6.2 At the end of the work day

 **NOTICE!** Do not leave any instruments on the motor overnight, in order to prevent oil from penetrating into the electric motor. Never lubricate the electric motor.

7 Care and maintenance

7.1 Cleaning and maintaining the mechanics using spray

Intervals


- At least daily at noon and in the evening
- Prior to sterilization
- After thermodisinfection


Required accessories

- | | |
|---|--------------------------|
| A | Spray adapter |
| B | O-ring for spray adapter |
| C | BA Ultimate Spray |

Procedure

- ✓ The motor has come to a stop.

✓ The O-rings on the spray adapter are intact.

1. Remove the file.
2. Pull the instrument off the motor.
3. Fit the spray adapter onto the nozzle of the spray can as shown.
4. Insert the instrument until it snaps into place and hold it firmly.
5. Apply spray to the instrument for 1 - 2 seconds.

Check

1. Check whether the liquid escaping from the instrument head is clean.
2. If the liquid is dirty: repeat the spraying procedure.

The liquid is still not clean?

1. Attach the instrument to the motor and let it run briefly. In this way, the spray can be distributed more effectively.
2. ⚠ CAUTION! Stop the motor.
Wipe off any leaking oil with a dry cloth.
3. Repeat the spraying procedure.

7.2 Cleaning the pushbutton chuck


To remove deposits and ensure proper functioning of the clamping system, the pushbutton chuck must be cleaned with BA Ultimate Spray.


Interval

Clean the pushbutton chuck at least once a week.

Procedure

1. Press the contra-angle handpiece head with the chuck firmly against the spray can nozzle.
2. Spray the chuck for 1 - 2 seconds.
Tip: Hold the spray can upright.
3. Wipe off any leaking protective oil with a lint-free cloth.


7.3 Manually cleaning and disinfecting the surface

✓ Wear appropriate protective clothing.

1. Remove the file [→ 17].
2. Remove the clip.
3. Automated processing is preferable to manual.

7.3.1 Cleaning the outer surface

⚠ NOTICE! Never clean in an ultrasonic bath!

1. Clean the instrument and the clips by brushing them off under running water (< 38°C, < 100°F, at least drinking water quality).
2. Blow it out with max. 3 bar.
3. Apply spray to the instrument [→ 24].

7.3.2 Disinfecting the surface

⚠ NOTICE! Never immerse in disinfectants!


✓ The disinfectants that are approved in your country must have proven bactericidal, fungicidal and virucidal properties.


1. Spray the surface with disinfectant.
2. Wipe the disinfectant away using a cloth.

The following can be used:

- MinutenSpray classic, by Alpro
- MinutenWipes, by Alpro

Please observe the manufacturer's instructions for using instrument disinfectants.

7.4 Sterilizing

1. Clean and disinfect the instrument [→ 35].
2. Apply spray to the instrument [→ 31].
3. Sterilize the instrument in the steam sterilizer with saturated water vapor.

Saturated water vapor: Temperature: 134°C (274°F)
Overpressure: 2.04 bar (29.59 psi)

Holding time: 3 min

Steam sterilizers that fulfill the requirements of either EN 13060, class B or EN 13060, class S and are suitable for the sterilization of straight/contra-angle handpieces are approved.

NOTICE! Do not exceed 140°C (284°F), even during the drying phase.

The instrument can be sterilized in a packaging which is suitable for sterilization and storage, e.g. a paper/plastic composite packaging.

After sterilizing

1. Remove the instrument from the steam sterilizer immediately.
 **CAUTION!** The instrument will be hot. There is a risk of burning!

NOTICE! Do not accelerate the cooling process by dipping the instrument into cold water. This will damage your instrument.

2. Store all instruments so that they are protected from contamination.
3. Sterilize again once the storage period has elapsed.

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